

tote-l-vac[®]

DESCRIPTION AND APPLICATION

Tote-L-Vac is a battery powered high vacuum and high flow suction unit. It features a soft, yet durable carrying case with easily accessible controls and battery. Although this unit is very compact, it provides high vacuum and high flow rates necessary for emergency oral aspiration procedures. Tote-L-Vac accepts both 800cc and 1200cc standard collection canisters. Tote-L-Vac is the perfect battery powered portable suction unit for EMS, home health, hospital patient transport and where ever battery power is a requirement.

FEATURES AND BENEFITS

- Vacuum regulator and gauge
- Housed in a durable carrying case
- Easy to clean and process for quick turnarounds
- Compact and lightweight
- **3 year warranty**

IDEAL FOR

It is an excellent choice for EMS, hospital patient transport, extended care facilities, doctors offices, home healthcare or where ever a battery powered suction unit is required.



BATTERY POWERED ASPIRATOR

E.M.S. • Patient Transport • Home Health Care • Surgicenters •
Nursing Homes • Physician's offices

TECHNICAL SPECIFICATIONS**Performance**

Vacuum Range: up to 550mmHg
 Free Air FLOW: 36 LPM
 Pump: 12 VDC oil-less diaphragm type

Collection Device

Canister: disposable plastic with mechanical shutoff
 Capacity: 800cc standard
 Tubing: 6' disposable sterile plastic
 Suction Tip: disposable Yankauer with control vent

Controls

Vacuum Regulator: rotary type on panel
 Vacuum Gauge: Calibrated to 0-30"/760 mmHg

Electrical Requirements

DC: 12 volts (nominal)

Battery:

Type: rechargeable sealed lead-acid
 Capacity: 12 volts, 3.0 ampere-hour
 Average Run Time: 30 minutes when fully charged
 Charger: fast/float charger
 UL/1310 listed
 Charge Time: 8 hours or less to 85% charging,
 depending on initial charge of
 battery
 Charge Method: 110 volt outlet, or 12 VDC vehicle
 outlet with optional cord

Physical

Weight: 15.5 lbs
 Dimensions: 8"H x 12"W x 9.5"D

PURCHASING SPECIFICATIONS

tote-I-vac® Part Number: 756000

SHIPPING INFORMATION

Shipping weight:
 20 lbs.

Package Size:
 12" H x 8" W x 8" D

Packaging:
 One per carton

Terms:
 net 30 days

Available Accessories

- Collection Canisters (disposable & glass)
- Aspirating Change Out Kits
- Filters & tubing
- Additional Accessories



Ohio Medical Corporation
 1111 Lakeside Drive
 Gurnee, IL 60031-4099 USA

www.ohiomedical.com

Form No. 550150 (Rev.7) 12 05

